A decorative graphic consisting of a trail of small, light blue particles that curves across the top half of the slide, resembling a DNA strand or a molecular path.

Revolutionizing Nucleic Acid Synthesis with Engineered Enzymes

David Entwistle Ph.D., Sr. Director, Codexis Inc.

TIDES EU November 1st 2023

CODEXIS[®]

ECO Synthesis™ Platform: Positioned to Deliver in RNAi Market

RNAi Demand is Coming




Chemical synthesis (phosphoramidite chemistry) alone **cannot meet projected future wave of demand for RNAi therapeutics**



Customers are asking us for a **scalable, more sustainable enzymatic solution** to complement chemical synthesis


1 kg of siRNA requires ~ 1000 kg of MeCN (BioSpring)

PMI per nucleotide added ~200 kg/kg (*Org. Proc. Dev*, 2021, 86, 1, 49–61)



Codexis positioned to deliver based on **20-year history of enzyme engineering** and **directly relevant Pharmaceutical Manufacturing commercial expertise**

Codexis ECO Synthesis™ Technology



Enzyme Catalyzed Oligonucleotide Synthesis

Core process:


Oligonucleotide synthesis by iterative, single nucleotide extension

- **Extend:** Controlled, enzymatic addition of modified, 3'-phosphate blocked ribonucleotides
- **Deblock:** Enzymatic cleavage of 3'-phosphate blocking group & excess ribonucleotides
- ... Repeat

Supply processes:

- Enzyme cascade for synthesis of modified, 3'-phosphate blocked ribonucleotides (NQPs)
- Enzymatic synthesis of starter oligonucleotide

Codexis ECO Synthesis™ Technology - Key Platform Traits



Enzyme Performance


- High incorporation efficiency
- Robustness & manufacturability

Scalable & Economical Manufacturing

- High volumetric productivity
 - “Oligo in solution, enzyme immobilized”
- Controlled addition of monomers
 - Low impurity production
- Smaller infrastructure/facilities footprint
- Established reagents supply

Supplying Critical NQP Reagents For ECO Synthesis™ Platform

Building a “one-pot, two-step” enzyme cascade with engineered kinases



Targeted Key Performance Indicators


Robustness/Manufacturability - soluble expression & stability

Productivity - >98% conversion at millimolar substrate concentration

Substrate Tolerance - accepts ribonucleotides with 2' modified sugars & phosphorothioate backbone modifications

Supplying Critical NQP Reagents For ECO Synthesis™ Platform

Building a “one-pot, two-step” enzyme cascade with engineered kinases




Step 1: N→NTP conversion via three consecutive phosphorylation steps, using three kinases

	A	C	G	U		A	C	G	U			
Initial (wild type) kinases	2'-OH	99	22	1	7	➔	2'-OH	98	96	93	91	Current engineered kinases
	2'-dF	64	0	0	0		2'-dF	93	19	93	29	
	2'-OMe	0	0	0	0		2'-OMe	55	0	0	0	

Percent conversion for individual nucleosides to the corresponding NTPs

Supplying Critical NQP Reagents For ECO Synthesis™ Platform

Building a “one-pot, two-step” enzyme cascade with engineered kinases




Step 1: N→NTP conversion via three consecutive phosphorylation steps, using three kinases

Status: NTP Forming Cascade

- Full base promiscuity & emerging activity for 2'-modified nucleosides
- Operational at process-relevant substrate concentration
- **Engineering on track to deliver full base and 2'-modification promiscuity for NTP formation**

Supplying Critical NQP Reagents For ECO Synthesis™ Platform

Building a “one-pot, two-step” enzyme cascade with engineered kinases



Step 2: NTP→NQP percent conversion with current 3'-O kinase


	A	C	G	U
Initial (wild type) kinase	2	0	0	0
Current engineered kinase	99	18	38	36

	A	C	G	U
2'-OH	99	18	38	36
2'-dF	12	0	0	0
2'-OMe	0	0	0	0

Percent conversion for individual NTPs to the corresponding NQPs

Supplying Critical NQP Reagents For ECO Synthesis™ Platform

Building a “one-pot, two-step” enzyme cascade with engineered kinases



Step 2: NTP → NQP percent conversion with current 3'-O kinase


Status: 3'-O-kinase engineering

- Engineering challenge as desired activity is non-natural (& potentially cytotoxic to expression host)
- Break through from initial strict A selectivity
- **Ongoing engineering of 3'-O-kinase starting to deliver activity on 2'-modification for NQP formation**

Terminal Deoxynucleotidyl Transferase (TdT)

A high engineered enzyme for catalyzing the extension reaction

PDB ID: 4I27



Targeted Key Performance Indicators

Robustness/Manufacturability - soluble expression & stability

Productivity - >99% conversion at millimolar substrate concentration

Substrate Tolerance - accepts ribonucleotides with 2'- and 3'-modifications & backbone phosphorothioate modifications


Promiscuity – minimal oligonucleotide sequence bias

Engineering Terminal Deoxynucleotidyl Transferase (TdT)

Progressive performance improvements in oligonucleotide synthesis

3'-Terminus Sequence	NQP	Iterative Rounds of Evolution							
		Starting TdT	1	2	3	4	5	6	7
AmC*mA*mG	mATP-3P	1	1	1	15	35	32	34	46
mAmAmAfG		2	2	2	2	2	3	4	58
mAmAfGmA		0	0	0	1	9	65	77	73
mAfGmUfG		0	0	0	0	0	0	0	39
...									
fGmUfGmU	mUTP-3P	0	0	0	1	31	48	67	68
mUfCmUfC		0	0	0	0	2	1	5	48
mUfCmAmU		0	0	0	3	55	35	56	78
fCmAmUmC		0	0	0	1	12	72	75	75
...									
mCmUmUmA	fATP-3P	0	0	0	2	10	76	65	68
mAmAmA(MOE)A		1	1	2	2	2	14	39	73
...									
mCmCmU*mU	mCTP-3P	0	0	0	2	9	47	66	68
	*mUTP-3P	0	0	0	0	1	33	54	69
	mGTP-3P	0	0	0	1	15	50	72	75
...									
*fA*fAfGmA	mCTP-3P	0	0	0	35	34	13	35	44
	*mUTP-3P	0	0	0	2	8	2	19	32
	mGTP-3P	0	0	0	35	41	20	42	52
...									



- One of the most complex evolution projects undertaken by Codexis to date
- Enabling synthesis of fully modified RNA oligonucleotides (2'-OMe, 2'-F, Phosphorothioate)



September 2023

Engineering Terminal Deoxynucleotidyl Transferase (TdT)



Progressive performance improvements in oligonucleotide synthesis



- Incorporation efficiency (N+1) of each NQP measured with 5-8 fully modified oligonucleotide substrates
- Included oligo modifications : 2'-OMe, 2'-F, Phosphorothioate)
- **Current process conditions:**
- [oligo] = mM; [NQP] = 2-fold eq
- **Ongoing enzyme engineering & process development**

Engineering Terminal Deoxynucleotidyl Transferase (TdT)

Progressive performance improvements in oligonucleotide synthesis




- Incorporation efficiency (N+1) on each oligonucleotide measured with 5 modified NQPs
- Included NQP modifications : 2'-MeO, 2'-F, Phosphorothioate)
- **Current process conditions:**
- [oligo] = mM; [NQP] = 2-fold eq
- **Ongoing enzyme engineering & process development**

Enzyme Performance: A Highly Engineered TdT

% Incorporation efficiency of N+1 additions over multiple rounds of evolution

- Consistent improvement in oligo promiscuity and acceptance of 2'-modified NQP
- Consistent movement of screening conditions to the desired process relevant targets
- Further engineered to remove remaining biases and increase oligo concentration to > 1 mM on going.
- Next steps, show process relevance

From Enzyme Engineering to Platform Development



Targeted Process Conditions

- Enzymes immobilized & oligos in solution
- Aqueous reaction system
- Substrates at millimolar concentration
- Elevated reaction temperature
- Cycle time per extension/deblocking: hours/NQP
- In-line monitoring of reaction progress


Platform performance

- Multiple consecutive cycles of oligonucleotide synthesis to produce a fully modified RNA sequence

Iterative oligonucleotide synthesis w/ ECO Synthesis™ platform

Platform performance

- 6+ extension cycles
- **Process conditions:**
oligo in solution, enzyme immobilized
- Extension with modified NQPs to yield fully modified RNA oligonucleotide
- Further enzyme engineering & process optimization ongoing



ECO Synthesis™ Technology: Controlled addition of monomers

- To our knowledge, the first ever multiple non-templated iterative extensions of an oligo in the solution phase mediated by highly engineered immobilized enzymes
- 6 extensions in batch mode with FAM-labelled oligo with immobilized enzymes
- Next step, demonstration in packed bed immobilized enzymes with solution phase unlabeled oligos in a flow set-up

ECO Synthesis™ Technology: A Vision For RNA synthesis

Unlocking exciting new possibilities in RNAi therapeutics manufacturing

- Accessing 10 - 100 kg batch size for oligonucleotide synthesis
- Reduction in organic burden & operational/facility costs
- Mild reaction conditions, elimination of protection groups, quality improvements

Codexis' deliverables

- Oligonucleotide synthesis process developed and optimized for product yield and quality
- Enzymatic manufacturing & supply chain for critical reagents (NQPs)
- Robust, high-performance enzymes for core & supply processes of ECO Synthesis™ platform, leveraging our CodeEvolver® enzyme engineering technology

Acknowledgements

ECO-Synthesis™ Technology - A company-wide cross functional collaboration

Behnaz Behrouzian

Margie Borra-Garske

Xingxiun Chen

Charlene Ching

Chinping Chng

Jackie Fleming

Stephanie Forget

Jeremiah Heredia

Callie Huitt-Roehl

Mikayla Krawczyk

Hirdesh Kumar

Michelle Li

Stefan Lutz

Melissa Mayo

Marissa MacAvoy

Niusha Mahmoodi

Matthew Miller

Michael Miller

Vesna Mitchell

Jovana Nazor

Simon Ng

Nick Porter

Philip Provencher

Asif Rahman

Gabriel Rajkovic

Ryan Reeves

James Riggins

Amani Shoubber

Leann Teadt

Ljubica Vojcic

Jonathan Vroom

David Watts

Leland Wong