
Improving Patients Subcutaneous Immunoglobulin (SCIg) Treatments Through Right Selection of Infusion Device Technology

Background

Previous studies have shown the importance of ancillary supplies on outcomes and corrective action dealing with local-site adverse events (AE's) during SClg administrations. The infusion device employed is a significant part contributing to the overall performance of all elements connected together in the delivery system to perform the infusion. There are three types of infusion devices available for SClg infusions. Each one creates different characteristics to support a successful SClg infusion outcome: electric type constant flow system (CFS), mechanical devices using variable coiled spring pressure (VSP) and mechanical constant pressure system (CPS).



Purpose

Identify the infusion device most appropriate for SCIg administration with considerations to the ancillary supplies, such as the subcutaneous needle set (SNS) and flow rate control device. Each system has advantages and disadvantages affecting the overall therapeutic outcome.

Contributing authors:


Fred L. Ma, M.D., Ph.D.¹, Chatarina Soderman, RN. MBA², Andrew Sealfon³, Paul Mark Baker M.D⁴

¹Chief Medical Officer, RMS Medical Products, ²VP of Product Marketing & Clinical Affairs Associate, RMS Medical Products,

³ President, RMS Medical Products, ⁴ Crystal Run Healthcare

Methods

Laboratory tests and clinical experiences have confirmed theoretical models of the three types of infusion devices and systems; CFS, VSP and CPS. The test in the laboratory included simulation of patient sites (depots) when different infusion devices where applied. Pressure at the outlet of the device or pump, and the distal pressure at the simulated patient depot were measured, recorded and compared.




Results

The results for the CFS show that as the depots saturate, the pressure at the simulated patient infusion-site and at the pump outlet increase. This process continues until the pressure (occlusion) alarm is activated.

The VSP shows decreasing pressure as fluid is delivered. This is a normal consequence of any coil spring as the force will drop with increasing distance. As a result, the measured pressure at the pump decreases during the delivery, resulting in decreasing flow rate.

The CPS shows a constant pressure at the pump independent of the volume delivered or the patient infusion-site pressure. The CPS reduced the flow rate in direct proportion to the simulated tissues saturation. For mechanical pumps and devices the flow rate is related to the differences in pressure between the pump and the patient depot.



Conclusion

The CPS demonstrated significantly lower pressures than CFS and more consistent flow rate performance and pressure control than VSP. For SCIg administrations, the CPS response to infusion-site saturation provides an improved level of safety and comfort, reduction of AE's and increasing patient compliance. AE's can be minimized and reduced through lower infusion pressures, defined and controlled flow rates limiting delivery when tissue perfusion is compromised as the patient's infusion-sites are saturating.

CPS is the system of choice:



- Significantly lower pressure than CFS
- More consistent flow rate performance and pressure control than VPS
- Responds to infusion-site saturation and provides feedback to improve level of safety and comfort
- Adverse effects can be reduced with:
- Lower infusion pressures
- Defined/controlled flow rates

References:

- 1. S. Jolles, J.S Orange, A. Gardulf, M.R. Stein, R. Shapiro, M.Borte and M.Berger, 2014, Current Treatment Options with Igg for the
- Individualization of Care in Patients with PID, British Society for Immunology, Clinical and Experimental Immunology 179: 146-160.

 2. Ochoa D, Curtis C, Duff C, Riley P, Murphy E and Zampelli A. 2012. Importance of Ancillary Supplies for Subcutaneous Immunoglobulin
- Infusions: Management of the Local Infusions Site, Journal of Clinical Immunology 32 (Supplement 1): 347.
- 3. Elena A. Skryabina, Teresa S Dunn, 2006, Supporting Documentation on Methods of Infusion, http://www.medscape.com/viewarticle/537694.
- 4. S. Keay and C. Callander, 2004, The Safe Use of Infusion Devices, Continuing Education in Anaesthesia, Critical care & Pain, Volume 4
- 5. Immunoglobulin National Society, IgNS Nursing Standard of Practice, First edition 2015.
- 6. B. Mainzer and K.D Stuhmeier, 1987, Aspekte des Druckaufbaus beim einsatz elektronischer infusionsgeräte, Anästh. Intensivher. Notfall med. 22, p.181 184.

TRS_PO2_InfusionDevices_IgNS2016_vA

