

Rohan T. Baker, Snigdha Gaur, David Murray, Aidan McEvoy, Graeme Young*, Lawrence LaPointe, and Susanne K. Pedersen Clinical Genomics Pty Ltd, North Ryde, NSW, Australia. *Flinders Centre for Innovation in Cancer, Adelaide, SA, Australia

BACKGROUND

We have previously evaluated the accuracy of a novel DNA blood test across a spectrum of benign and neoplastic conditions in the colon/rectum using a multiplexed methylation-specific PCR assay for detection of circulating <u>fully</u>-methylated BCAT1 and IKZF1 DNA (Assay 1). However, a number of methylation patterns are likely to exist in a colorectal neoplasm due to the heterogeneous nature of human solid tumours.

AIM

To determine if a modified 2-gene blood test that enables detection of partial methylation of three interprimer CpG sites in the IKZF1 target region (Assay 2) improves detection of cancer.

MATERIALS and **METHODS**

Cell-free DNA recovered from 4mL plasma from colonoscopy confirmed subjects was bisulphite-converted and assayed for the presence of methylated BCAT1 and IKZF1 DNA (Pedersen et al., 2015). The BCAT1 MethyLight assay was multiplexed with one of two different *IKZF1* MethyLight assays (Assay 1 or Assay 2; see Assay Details below). The *IKZF1* PCR component in Assay 1 constituted three methylation-specific oligonucleotides (forward and reverse primers, plus a hydrolysis probe) for amplification and detection of a fully methylated 95-bp target region. The *IKZF1* PCR component in Assay 2 used the same two methylation-specific *IKZF1* primers as Assay 1, but used a 'degenerate' hydrolysis probe cocktail, synthesised with a 50:50 mix of C and T at three interprimer CpG sites. Thus, Assay 2 would detect eight different interprimer methylation patterns in the 95-bp *IKZF1* target region.

Assay Details

- BCAT1 Fwd primer: 5'-GTTTTTTGTTGATGTAATTCGTTAGGTC
- BCAT1 Rev primer: 5'- CAATACCCGAAACGACGACG
- BCAT1 Probe: HEX-5'-TTCGTCGCGAGAGGGTCGGTT-lbFQ
- IKZF1 Fwd primer: 5'-GACGACGTATTTTTTCGTGTTTC
- IKZF1 Rev primer: 5'-GCGCACCTCTCGACCG
- IKZF1 Probe (Assay 1): FAM-5'-TTTGTATCGGAGTAGCGATTCGGGAG-IbFQ
- IKZF1 Probe (Assay 2): FAM-5'-TTTGTATYGGAGTAGYGATTYGGGAGG-IbFQ (Y = 50% C, 50% T). • ACTB Fwd primer: 5'-GGAGTTTTTGTTTTTGGTTAGTTG
- ACTB Rev primer: 5'-CAAAATAAAATACAAAACAAACCTAATCC
- Tex615-5'-ATGGAGGTTTAGTGGTAATATAGGTTTTGTTTGG- IbRQ ACTB Probe:
- PCR conditions:
- 30µL final PCR volume; 12µL bisulfite-converted DNA; 15µL 2X Quantitect NOROX (QIAGEN); 3µL 10X oligo mix. All primers at 200nM final; all probes at 100nM final
- 1 cycle, 15min/95°C; 50 cycles: 15sec/95°C; 40sec/62°C with acquisition (FAM, HEX and Texas Red).
- Roche LC480 II with colour-compensation applied.

Detection of variable methylation patterns improves sensitivity of a colorectal cancer blood test

ASSAY 1 POSITIVITY

Assay 1 positivity rate by clinical status

Clinical status	n	CVP 2-gene panel		
Cinnear Status	2108	n	%	
Non-neoplastic	1283	74	5.8	
Non-advanced Adenomas	460	30	6.5	
Advanced Adenomas	232	17	7.3	
TIs (stage 0)	3	0	-	
Cancers	130	85	65.4	
Stage I	24	7	29.2	
Stage II	53	36	67.9	
Stage III	39	30	76.9	
Stage IV	9	8	88.9	
Unstaged	5	4	80.0	

65% sensitivity (any cancer) / 94% specificity



- *IKZF1* amplicons were sequenced from bisulphiteconverted cfDNA extracted from colonoscopy-confirmed cancer cases that were positive (left panel) or negative (centre and right panel) with the fully-meth *IKZF1* probe.
- The centre and right panels give examples of cancers that show variable methylation at the *IKZF1* locus.
- The *IKZF1* assay was thus redesigned to include a mixture of 8 probes to detect all possible methylation combinations (See IKZF1 Probe (Assay 2), Assay details).

References

- Pedersen SK, Baker RT, McEvoy A, Murray DH, Thomas M, Molloy PL, Mitchell S, Lockett T, Young GP, LaPointe LC. A Two-Gene Blood Test for Methylated DNA Sensitive for Colorectal Cancer. PLoS One. 2015 Apr 30;10(4):e0125041. doi: 10.1371/journal.pone.0125041.
- Young, GP, Pedersen, SK, Dekker, E, Cole, SR, Osborne, JM, Symonds, EL, Mallant-Hent, RC, Mcevoy, A, Baker, R, Gaur, S, Murray, DH, Lapointe, LC. Evaluation of a 2-gene (IKZF1 and BCAT1) DNA blood test for detection of colorectal cancer. Digestive Disease Week, Chicago 2014. Gastroenterology 2014;79(5) Suppl AB125

Young et al., DDW 2014

Rev primer Partially methylated

RESULT 2: ASSAY 2 POSITIVITY Assay 2: *IKZF1* positivity rate

Clinical status	Ν	Full-meth		Vari-meth		McNemar
	743	n	%	n	%	Р
Non-neoplastic	514	8	1.6	17	3.3	0.08
Adenomas	196	3	1.5	11	5.6	0.04
TIs (stage 0)	2	0	-	1	50	N/A
Cancers	33	11	33.3	17	51.5	0.04
Stage I+II	24	6	25	11	45.8	0.07
Stage III+IV	9	5	55.6	6	66.7	1.00

- A subset of 743 previously assayed samples (See panel at left) were re-assayed with Assay 2, allowing detection of variably-methylated IKZF1.
- *IKZF1* positivity increased for all cancers from 33.3% to 51.5% (p=0.04), and particularly for early stage (I+II) cancers from 25% to 46% (p=0.07).

RESULT 3: ASSAY 2 CLINICAL PERFORMANCE

Assay 2 in a clinical cohort

2014: Partial IKZF1 meth assay 2013: Full-meth assay n = 677 (467N, 175A, 2 Stage 0, 33 CRC)



64% sensitivity/94% specificity

• When combined with BCAT1 results, Assay 2 resulted in an increase from 64% to 70% sensitivity for cancers, with a 92% specificity.

 The increased sensitivity was due to an increase in positivity in early stage cancers.

CONCLUSIONS

- Allowing detection of variably-methylated *IKZF1* improved blood-test sensitivity for early stage cancers.
- *IKZF1* methylation may be incomplete in early cancers.
- Seeking opportunities to further investigate clinical utility.
 - info@clinicalgenomics.com

ColoVantage[®]

Live life well.



70% / 92%