

# MEDICAL GRADE SILICON NITRIDE



A M E D I C A<sup>®</sup>

Ametica's proprietary composition of silicon nitride provides the right combination of strength, toughness, wear resistance, biocompatibility, bioactivity, bone integration, structural stability, corrosion resistance, and easier imaging, all of which are desirable in medical implants.



## THE IDEAL BIOMATERIAL

**Strength and fracture toughness:** Interlocking anisotropic grains deflect and bridge cracks.

**Wear resistance:** High hardness, strength, and fracture toughness prevent wear.

**Material phase stability:** No spontaneous phase transformation or associated weakening.

**Hydrophilicity:** Tunable through modification of surface topography and chemistry, from  $<10^\circ$  up to  $\sim 70^\circ$ .

**Osseointegration:** Nanostructured topography combined with complex surface chemistry optimal for cell adhesion and bone growth.

**Favorable imaging:** Semi-radiolucent density appears bone-like in X-rays and low magnetic susceptibility eliminates distortion in CT and MRI scans.

**Bacterial resistance:** Surface chemistry, nanotexture, and surface charge inhibit biofilm formation.

## TYPICAL PROPERTIES OF AMEDICA'S SILICON NITRIDE

Property	Units	Test Method	Typical	Specification
Density	g/cc	ASTM C 373	3.26	≥ 3.23
Grain Size	microns	BS EN 623-3	0.5 x 5.0	≤ 25
Flexural Strength	MPa	ASTM C 1161	1,000	≥ 900
Weibull Modulus	-	ASTM C 1239	10	≥ 6
Elastic Modulus	GPa	ASTM C 1161	296	≥ 290
Poisson's Ratio	-	*	0.27*	-
Compressive Strength	MPa	*	>4,000*	-
Hardness	GPa	ASTM C 1327	15.0	≥ 14.3
Fracture Toughness	MPa·m <sup>1/2</sup>	ASTM E 399	10.5	≥ 9.0
Si3N4 Phase Composition	%	Xray Diffraction JCPDS# 82-0697	100% β-Si3N4	≥95% β-Si3N4
Biocompatibility	-	ISO 10993	Pass	Pass
Thermal Conductivity	W/m·°K	*	15-30*	-
Coefficient of Thermal Expansion (RT – 200°C)	1 x 10 <sup>-6</sup> /°C	*	2.2*	-
Specific Heat	J/Kg·°K	*	170*	-
Volume Resistivity	ohm-cm	*	>1012*	-

\*Reported data are typical of silicon nitride. These values have not been specifically measured for Amedica's silicon nitride.

## BIOCOMPATIBILITY TESTING

Test	Method
Cytotoxicity	ISO 10993-05
Sensitization	ISO 10993-10
Intracutaneous Toxicity	ISO 10993-10
Acute Systemic Toxicity	ISO 10993-11
Subchronic Toxicity	ISO 10993-11
Genotoxicity	ISO 10993-3
Muscle Implant tests	ISO 10993-6
Physicochemical Testing	USP

## COMPATIBILITY WITH STERILIZATION METHODS

	Gamma Irradiation, E-Beam, X-Ray	Steam	Ethylene Oxide Gas
Silicon Nitride	Yes	Yes	Yes

## ASTM OR ISO SPECIFICATIONS FOR BIOMATERIALS TABLE<sup>1</sup>

Property	Si <sub>3</sub> N <sub>4</sub> ASTM F-2094 <sup>1</sup>	Si <sub>3</sub> N <sub>4</sub> ISO 266021	Al <sub>2</sub> O <sub>3</sub> ASTM F-603	Al <sub>2</sub> O <sub>3</sub> ISO 6474-1	Mg-PSZ ASTM F-2393	Y-TZP ASTM F-1873	ZTA, AMC ISO 6474-2	CoCr ASTM F799	Ti6Al4V ASTM F136	PEEK ASTM F2026
Chemical Purity (%)	≥ 97.0	NS	≥ 99.5	≥ 99.7	≥ 99.8	≥ 99.0	≥ 99.8	NA	≥ 99.3	NA
Density (g/cc and %)	3.0 – 3.4 ≥ 99.8	3.0 – 3.6 NS	≥ 3.93 ≥ 98.6	≥ 3.94 ≥ 98.8	≥ 5.80 ≥ 98.8	≥ 6.00 ≥ 98.4	≥ 4.31 ≥ 98.6	NA	NA	1.28 – 1.32
Grain Size (µm)	NS	NS	≤ 4.5	≤ 2.5	NS	≤ 0.6	Al <sub>2</sub> O <sub>3</sub> ≤ 1.5 ZrO <sub>2</sub> ≤ 0.6	≤ 64	NA	NA
Flexural Strength (MPa) <sup>1</sup>	≥ 765	≥ 760	≥ 400	≥ 500	≥ 600	≥ 800	≥ 750	827 (YS) 1172 (TS) <sup>2</sup>	760 (YS) 825 (TS)	110
Weibull Modulus	≥ 12	≥ 12	≥ 8	≥ 8	≥ 10	NR	≥ 8	NA	NA	NA
Fracture Toughness (MPa·m <sup>1/2</sup> )	≥ 6.0	≥ 6.0	NS	≥ 2.5	NS	NS	≥ 3.5	NA	NA	50 <sup>3</sup>
Hardness (GPa)	≥ 15	≥ 14.2	≥ 18	≥ 18	≥ 10	≥ 12	≥ 15.5	≥ 3.3	NA	NA
Elastic Modulus (GPa)	270 – 330	270 – 330	≥ 380	≥ 380	≥ 180	≥ 200	≥ 320	NA	NA	3


## PHYSICAL & MECHANICAL PROPERTIES AND PERFORMANCE OF BIOMATERIALS<sup>1</sup>

Property or Performance	Units	Silicon Nitride	Alumina	Zirconia	Zirconia-Alumina Composites	Cobalt Chromium	Ti6Al4V	PEEK	Cortical Bone	Titanium Nitride	Diamond-Like Carbon	Zirconium Nitride	Titanium Niobium Nitride	Oxidized Zirconium	Hydroxyapatite	
Composition	NA	Si <sub>3</sub> N <sub>4</sub>	Al <sub>2</sub> O <sub>3</sub>	Mg-PSZ Ce- or Y-TZP	m-ZTA AMC ATZ	ASTM F799	ASTM F136	ASTM F2026	Collagen, Proteins, HAp	TiN	DLC	ZrN	TiNbN	Ox-Zr	ASTM F1609	
Density	g/cc	3.22 – 3.35	3.98	5.65 – 5.77 6.00 – 6.05	4.25 4.37 5.51	8.29 – 8.50	4.43 – 4.50	1.29	1.5 – 2.0	4.87 – 5.22	0.90 – 3.20	7.09	~5.69	5.84	2.55 – 3.21	
Grain Size	µm	0.5 x 5.0 Non-Equiaxed	<1.8 Equiaxed	50 Equiaxed 0.1 – 0.6 Equiaxed	0.4 – 0.7 Equiaxed 0.54 Equiaxed 0.4 Equiaxed	~62 Equiaxed	~10 x 60 Lamellar	NA	NA	30 – 300 nm Columnar	Amorphous 2 - 25 nm	10 – 30 nm Nanocrystals	10 – 30 nm Nanocrystals	40 x 200 nm	0.4 - 100 µm splats	
Flexural or Tensile Strength	MPa	800 – 1100 Flexural	400 – 580 Flexural	450 – 700 Flexural 700 – 1500 Flexural	700 – 1248 Flexural 1250 – 1400 Flexural 755 – 1163 Flex./Biaxial	827 Tensile	860 – 970 Tensile	170 Flexural	90 – 228 Flexural	10 – 60 N IC Adhesion	35 – 160 N IC Adhesion	24 – 60 N IC Adhesion	83 N IC Adhesion	35 N IC Adhesion	39 – 189 Bond 25 – 60	
Compressive Strength	MPa	4000	4100 – 5000	2000 - 3000 2000 – 2200	4000 - 4500 4300 ~2600	600 – 1800	800 - 970	118	150 – 260 // 70 – 110 ⊥	400 – 5500	NA	NA	NA	~2000	102 – 1000	
Elastic Modulus	GPa	296 – 313	380	200 – 250 210 – 223	340 – 390 358 240 - 250	197 – 210	105 – 120	4	7.5 – 25.8 // 5 – 20 ⊥	402 – 550	110 – 900	175 – 395	200 – 600	200	3.2 to 122 coat vs. bulk	
Poisson's Ratio	NA	0.27	0.23	0.30 0.30 – 0.33	~0.24 0.24 ~0.28	0.27 – 0.32	0.31 – 0.34	0.4	0.19 – 0.48	0.21	0.17 – 0.20	0.19	~0.20	0.34	0.11 – 0.27	
Weibull Modulus	NA	8 – 53	5 – 29	22 7 – 87	NA 10 – 15 6 – 17	NA	NA	NA	NA	5 – 18	6 – 12	NA	NA	NA	2 – 19	
Fracture Toughness	KIC, MPa·m <sup>1/2</sup>	4.4 – 15.0	3.3 – 4.2	2.9 – 16.0 4.5 – 20.0	>4.1 6.4 – 8.5 8.0 – 12.0	50 – 100	46.3 – 93.3	7.6 kJ/m <sup>2</sup> Impact Test	1.0 – 5.0 // 3.0 – 20.0 ⊥	0.7 – 12.4	1.6 – 5.1	2.3 – 7.5	NA	2.2 – 2.8	0.5 – 1.2	
Fatigue Resistance	K <sub>TH</sub> /K <sub>C</sub>	0.50 – 0.97	0.52 – 0.84	0.45 – 0.90 0.37 – 0.92	NA 0.67 NA	0.14 – 0.36	0.10 – 0.40	0.53 – 0.62	0.30 – 0.83	NA	NA	NA	NA	NA	0.61	
Biocompatibility	NA	Pass	Pass	Pass Pass	Pass Pass Pass	Marginal	Pass	Marginal	Pass	Pass	Pass	Pass	Pass	Pass	Pass	
Surface Phase Composition	%	100% β-Si <sub>3</sub> N <sub>4</sub>	100% α-Al <sub>2</sub> O <sub>3</sub>	42% - 54% t-ZrO <sub>2</sub> 65% - 95% t-ZrO <sub>2</sub>	83% - 93% t-ZrO <sub>2</sub> 58% - 90% t-ZrO <sub>2</sub> 95% - 99% t-ZrO <sub>2</sub>	NA	Mixture of α & β Ti	Amorphous & Crystalline	Collagen and HAp	TiN Nanocrystals	Amorphous	ZrN Nanocrystals	TiN, NbN Nanocrystals	95% m-ZrO <sub>2</sub> 5% t-ZrO <sub>2</sub>	ACP, TCP, HA	
ITD Susceptibility	NA	Stable	Stable	Marginal Metastable (Y-TZP; Marginal (Ce-TZP)	Stable Marginal Metastable	Stable	Stable	Stable	Stable	Stable	Stable	Stable	Stable	Stable	Purposely Unstable	
Hardness	GPa	15.0	18.0 – 23.0	10.0 – 12.0 11.0 – 12.5	15.7 – 20.8 19 13.7 – 15.0	3.0 - 4.0	2.8 – 3.3	99 Rockwell M	0.68 – 0.78 // 0.46 – 0.57 ⊥	33 – 56	14.5 – 80.0	14.0 – 31.0	14.0 - 24.5	12.0 – 14.0	3.0 – 9.0	
Wear Rate PE HXLPE Hard-on-Hard	mm <sup>3</sup> /MC	17 - 25 3.7 – 6.3 0.18 – 0.98	20 – 58 0.0 – 6.9 0.02 – 4.71	- 1.8 – 5.1 NA	11 – 63 5.0 – 6.0 Catastrophic	NA	1 – 20 0.1 – 4.4 0.00 – 0.45	17 - 32 5.6 – 6.1 0.02 – 0.06	14 – 201 0.0 – 11.7 0.18 – 25.00	NA	NA	NA	21 NA NA	28 – 67 2.8 NA	NA 3.5 NA	NA
Thermal Expansion Coefficient	10 <sup>-6</sup> /°C	2.0 – 4.6	8	7 - 10 11	~8 8.1 ~10	7.32	8.5 – 9.7	47	22.0 – 32.4	7.4 – 9.2	2.3	5.9 - 7.2	~7.4 – 9.2	7 – 10	11.6 – 14.2	
Thermal Conductivity	W/m·°K	30 - 50	30	2 2 – 3	~17 17 ~6	12.7	6.7 - 7.0	0.29	0.41 – 0.63	11.9	0.2 – 30	20	~12 - 14	2 – 3	1.1 – 1.2	
X-Ray Radiolucency	NA	Radiolucent	Radiolucent	Opaque Opaque	Opaque Opaque Opaque	Opaque	Opaque	Transparent	Radiolucent	Opaque	Opaque	Opaque	Opaque	Opaque	Radiolucent	
Sessile Water Contact Angle	Degree (°)	40 – 70	50 – 72	79 82	90 90 90	55 – 93	76	95	NA	31 - 69	55 – 71	89	73 - 75	71	34 – 39	
Bacteriostatic Capabilities	⊕ = Excellent; + = Good; ∅ = Fair; ⊗ = Poor; x = Very Poor	⊕	∅	+	+	NA	+	NA	⊗	+	x	NA	+	+	+	
Osseointegration Ability		⊕	∅	+	+	NA	+	NA	⊗	+	x	⊕	+	+	⊕	

<sup>1</sup> B.J. McEntire, B.S. Bal, M.N. Rahman, J. Chevalier, and G. Pezzotti, "Ceramics and Ceramic Coatings in Orthopaedics," J. Eur. Ceram. Soc., **35** [16] 4327-4369 (2015).


## SILICON NITRIDE FORMS

### As Fired




Nano-scale topography, increased surface area  
**Ideal for:** Bone on growth

### Textured




Customizable macro-scale topography  
**Ideal for:** Grip, bone on growth

### Porous




Up to 70% connective porosity with 100-700  $\mu\text{m}$  pore size  
**Ideal for:** Bone ingrowth

### Polished



$R_a \sim 5 \text{ nm}$ , Hardness  $\sim 15 \text{ GPa}$   
**Ideal for:** Bearing surfaces

### Composite



Application-specific custom morphologies  
**Ideal for:** Hybrid applications

## SHARING OUR EXPERTISE

At Ametica, we have the scientific and manufacturing expertise to produce medical grade silicon nitride - a patented platform technology for use in a variety of medical applications. Silicon nitride is bioactive and compatible across all imaging modalities, offering surgeons and patients a preferable alternative to commonly used materials.

**Let our leading R&D and manufacturing teams convert your existing medical devices into silicon nitride. With our unrivaled in-house capabilities, we are equipped to control complex geometries on a macro-, micro-, and nano-level, which allows for intricate designs and shapes that can be rapidly developed, prototyped, and tested in our FDA registered and ISO 13485 certified facility.**

Contact us at **855.839.3500** or online at **amedica.com** to learn more.



© 2016 AMEDICA Corporation. All rights reserved.  
1885 West 2100 South, Salt Lake City, UT 84119  
855.839.3500 (toll free) | [www.amedica.com](http://www.amedica.com)

CALL **855.839.3500** OR VISIT US AT **AMEDICA.COM** TO DISCOVER WHY MATERIAL MATTERS.