

The background features several thick, curved green lines that sweep across the slide, creating a dynamic, organic feel. These lines are in various shades of green, from light to dark.

Engineered Enzymes to Overcome Scalability and Sustainability Challenges of Nucleic Acid Therapeutics Manufacturing

Mathew Miller, PhD Associate Director, Life Sciences Technology and Applications


Stephanie Forget, PhD Senior Scientist I

Codexis core business pillars

Based on CodeEvolver® platform to accelerate enzyme discovery and commercialization



Enzymes
from Nature



Commercially
Relevant Enzyme

Value
Creating
Products

Pharmaceutical Manufacturing
enzymes for small molecule
production

Life Sciences
enzymes for NGS applications
and DNA/RNA synthesis

Biotherapeutics
enzymes as oral drugs;
engineered transgenes and
capsids for gene therapy

About the Speaker




Mathew Miller, Ph.D.

Associate Director of Life Science Technology and Applications, Codexis

Mathew joined Codexis in 2009. He is currently focused on enzymes used in mRNA manufacturing, including directed evolution for T7 RNA polymerase. Matt has previously lead protein engineering projects for next generation sequencing and enzymatic DNA synthesis. He received his Ph.D. from the University of California, San Francisco.

mRNA is a highly versatile therapeutic modality with manufacturing challenges



Applications


- Vaccines (SARS-CoV2, others)
- Cancer vaccines (Immuno-oncology)
- Gene editing (CRISPR/Cas, others)
- Therapeutic proteins

mRNA manufacturing challenges

- mRNA capping efficiency
- Undesired side products (dsRNA)
- mRNA yield
- Downstream processing

Eukaryotic mRNA cap: structure and function

7-Methylguanosine (7mG) Cap:



Cap functions:


- Prevents recognition of mRNA as foreign/non-self
- Stabilizes mRNA
- Initiates translation

Co-transcriptional capping

Dinucleotide cap analog:

N7-Methyguanosine

Guanosine



T7 RNA Polymerase


competition
vs GTP



Cap analog is the most expensive reaction component

Increased ratio of cap to GTP ensures:

- Increased capping efficiency
- Decreased mRNA yield
- High cost of production



RNA Polymerase Evolution Target

Enzyme Engineering targets:

Increased capping efficiency




- Reduce cost of the process

Reduced dsRNA generation




- Increase mRNA safety and efficacy
- Reduce cost of the process

CodeEvolver® technology rapidly generates enzymes with desired properties through iterative rounds of enzyme optimization



Co-transcriptional capping: ARCA cap analog



ARCA




- Inefficiently incorporated by WT T7RNAP

Codex[®] HiCap RNA Polymerase

Increased dinucleotide analog capping efficiency (ARCA)



Optimized condition:




- Allows higher capping efficiency for a given concentration of ARCA, GTP
- Allows for >95% capping efficiency in batch IVTs without GTP starvation, preserving yield and lowering cost

Co-transcriptional capping: CleanCap® AG cap analog

CleanCap®
Reagent AG

Trinucleotide




- Higher incorporation by WT T7RNAP
- Native Cap-1 structure
- More expensive to produce and use


Codex® HiCap RNA Polymerase

Increased trinucleotide (CleanCap® AG) capping efficiency


Capping efficiency



Capping Efficiency




mRNA yield



- Codex® HiCap RNA Polymerase selectively incorporates CleanCap® AG
- Allows for higher capping efficiency or equivalent efficiency with lower cap concentrations (1.5mM), reducing cost

Double-stranded RNA side products






- dsRNA is generated by WT T7 via an RNA-dependent RNA polymerase activity
- Pattern recognition receptors (TLR3, RIG-I) recognize dsRNA and induce interferon signaling, causing inflammation
- dsRNA signaling inhibits protein translation
- dsRNA contamination impacts safety and efficacy of mRNA therapeutics
- Polishing step required to remove dsRNA – adds complexity and cost


Codex[®] HiCap RNA Polymerase

Reduced dsRNA formation

Analytical reverse-phase HPLC:




dsRNA-specific Sandwich ELISA:




- 3' extensions for the HiCap RNA Polymerase sample are significantly reduced vs the WT sample
- Residual dsRNA from HiCap RNA Polymerase without purification is similar to that from purified WT polymerase sample
- dsRNA sandwich ELISA shows ~50-fold reduction in total dsRNA side products
- Lower dsRNA eliminates the need for polishing – simplified workflow and lower cost

Reduced dsRNA formation correlates with reduced immune stimulation and increased protein translation

Reduced immune stimulation



GFP reporter translation post-transfection



- Codex® HiCap RNA Polymerase shows significant reduction in IFN signaling compared to WT T7
- Codex® HiCap RNA Polymerase shows significant increase in protein translation in HeLa cells
- Lower dsRNA improves safety and efficacy, using a simplified workflow

Codex[®] HiCap RNA Polymerase attributes

	Codex [®] HiCap RNA Polymerase	Wild-type T7 polymerase
Reduced dsRNA production	● ● ●	●
Co-transcriptional capping efficiency	● ● ●	● ●
IVT yields	● ● ●	● ● ●
Reduced capping reagent	● ● ●	●
Manufacturing simplification	● ● ●	●
Finished product yield	● ● ●	●

Information previously presented by Tyler Goodwin, Ph.D., Head of Non-viral Delivery Development and Vector Core, Tune Therapeutics

Higher quality mRNA and improved manufacturing cost & efficiency

The background features several thick, curved green lines that sweep across the slide, creating a dynamic, organic feel. These lines vary in shade from a light, pale green to a vibrant, medium green.

Engineered Enzymes to Overcome Scalability and Sustainability Challenges of Nucleic Acid Therapeutics Manufacturing

Mathew Miller, PhD Associate Director, Life Sciences Technology and Applications

Stephanie Forget, PhD Senior Scientist I

RNAi: High Demand / Constrained Supply

Demand Drivers



RNA delivering on promise of personalized medicines
(10+ FDA approved therapies have reached market in past 5yrs)

RNAi therapeutics as a modality is growing rapidly with
>450 assets in pipelines

RNAi therapies are treating the previously untreatable diseases

Production & Supply Challenges



Chemical RNA synthesis produces **Millions** of liters of
chemical waste... and growing!

Phosphoramidite synthesis will be challenged to meet
demand of 1,000s kg of RNAi p.a. by 2030

Critical solvent supply constraints likely (Acetonitrile)

Enzymes are poised to spur innovation and disrupt RNAi manufacturing

Scalability



Enzyme catalytic activity
enables multiple cycle use

Purity & Yield



Evolution targets >99%
incorporation efficiency

Sustainability



Aqueous waste streams &
lower solvent use

About the Speaker



Stephanie Forget, PhD

Senior Scientist I, Codexis


Technical lead for enzymatic RNA synthesis program. Obtained PhD in Chemistry at Dalhousie University, Canada. Studied biosynthetic enzymes involved in rare sugars found in natural products. Postdoctoral fellow at University of British Columbia, Canada. Studied oxidative enzymes through reaction profiling.



CODEXIS[®]


ECO Synthesis[™] Technology

(Enzyme-Catalyzed **O**ligonucleotide Synthesis)



Codexis ECO Synthesis™ Technology

Enzyme-Catalyzed Oligonucleotide Synthesis for RNAi therapeutics




ECO Synthesis™ Technology

- Controlled addition of modified RNA bases (TdT)
- Deblocking of 3' blocking group (phosphatase)
- Supply of 3' blocked NTP substrates (multiple enzymes)

Codexis ECO Synthesis™ Technology

Enzyme-Catalyzed Oligonucleotide Synthesis for RNAi therapeutics



Final process in development

Enzyme Performance

- High incorporation efficiency (>99%)
- No sequence bias

At-Scale Process Requirements


- Controlled addition of monomers
- Low impurity production
- High volumetric productivity

Scalable & Economical Enzyme Manufacturing

- High manufacturing yields
- Scalable supply of nucleotide triphosphates

Codexis ECO Synthesis™ Technology

Enzyme-Catalyzed Oligonucleotide Synthesis for RNAi therapeutics



Final process in development

Enzyme Performance

- High incorporation efficiency
- No sequence bias

At-Scale Process Requirements


- Controlled addition of monomers
- Low impurity production
- High volumetric productivity

Scalable & Economical Enzyme Manufacturing

- High manufacturing yields
- Scalable supply of nucleotide triphosphates

Codexis ECO Synthesis™ Technology

Enzyme-Catalyzed Oligonucleotide Synthesis for RNAi therapeutics



Final process in development

Enzyme Performance

- High incorporation efficiency
- No sequence bias

At-Scale Process Requirements


- Controlled addition of monomers
- Low impurity production
- High volumetric productivity

Scalable & Economical Enzyme Manufacturing

- High manufacturing yields
- Scalable supply of nucleotide triphosphates

Codexis ECO Synthesis™ Technology

Enzyme-Catalyzed Oligonucleotide Synthesis for RNAi therapeutics



Final process in development

Enzyme Performance

- High incorporation efficiency
- No sequence bias

At-Scale Process Requirements

- Controlled addition of monomers
- Low impurity production
- High volumetric productivity


Scalable & Economical Enzyme Manufacturing

- High manufacturing yields
- Scalable supply of nucleotide triphosphates

Enzyme Performance: A highly engineered TdT

Key Enzyme Challenges

- Does not naturally recognize modified RNA
- Poor soluble expression and stability
- Enzymes function close to physiological conditions
- Manufacturability for large scale enzyme needs



Enzyme Performance: A highly engineered TdT

% Incorporation efficiency of N+1 additions over multiple rounds of evolution

3'-Terminus Sequence	NTP	Starting TdT Iterative Rounds of Evolution								
mAmCmU	fATP-3P	0	0	0	0	1	54	91	94	
	fUTP-3P	0	0	0	0	1	56	93	94	
	mCTP-3P	0	0	0	0	1	54	91	89	
	mGTP-3P	0	0	0	0	1	24	87	78	
	mATP-3P	0	0	0	0	1	39	88	66	
	*mGTP-3P	0	0	0	0	0	0	2	55	
	*mUTP-3P	0	0	0	0	0	0	19	24	
mGmAmC	fUTP-3P	0	0	0	0	2	33	74	92	
	mATP-3P	0	0	0	0	1	16	60	90	
	fATP-3P	0	0	0	1	2	34	77	88	
	mCTP-3P	0	0	0	0	0	13	66	86	
	*mUTP-3P	0	0	0	0	0	0	11	79	
	*mGTP-3P	0	0	0	0	0	0	3	64	
	mGTP-3P	0	0	0	0	0	0	2	55	
AT*mG	mATP-3P	4	1	4	49	75	82	47	56	
AmU*mG		0	0	0	12	46	77	75	41	
mAmU*mG		0	0	0	0	5	69	82	68	
mAmUfG		0	0	0	0	0	2	1	66	
mUmGmA	mATP-3P	0	0	0	1	2	38	82	86	
mAfUCmC		0	0	0	0	4	58	88	86	
mAmG(MOE)C		0	0	0	0	0	4	14	84	
mC*mG*mA		0	0	0	0	2	39	75	75	
mCmUmG		0	0	0	0	4	82	86	72	
mAmUmC		0	0	0	0	2	57	84	63	
mAmUfU		0	0	0	0	0	0	0	59	
*mA*mA*mC		0	0	0	0	1	30	54	57	
mAmUfC		0	0	0	0	0	1	1	51	
mCmGmA	fATP-3P	0	0	0	1	4	82	93	92	
*fAfGmA		4	2	4	19	50	79	85	85	
mC*mG*mA		0	0	0	1	3	47	81	82	
*fA*fAfG		0	0	0	0	0	44	56	70	
fCfGfA		0	0	0	3	26	50	55	65	
mU*fA*fA		0	0	0	0	3	11	17	34	
fGmAfU		6	5	6	3	12	42	52	30	
fC*fG*fA		0	0	0	0	4	14	30	14	

So far...

- ✓ Incorporation of relevant 2'-modifications with 3' phosphate blocking group
 - ✓ 2'-deoxyfluoro and 2'-methoxy
- ✓ Incorporation of α-PS bonds
- ✓ Recognition of modified initiator sequences
- ✓ Improving incorporation efficiency of each nucleotide (goal >99%)

Note: All incorporating nucleotides contain 3'-blocking group; "*" denotes alpha PS-bond; "m, r, f" denotes 2'-OMe, 2'-OH, or 2'-F modifications, respectively




5' – Biotin/FAM...mA mU mA **mA** mG **mA** fA – 3'

- ✓ 4 Cycles for proof of concept
- ✓ Feasibility demonstrated for modified RNA synthesis

ECO Synthesis™ Technology: Controlled addition of monomers

Addition and deblocking with 5'-biotynilated-FAM...mAmUmA




Not final process conditions

Coming soon...enzymatically-activated monomer supply

Demonstrated...

✓ Nucleoside → NTP conversion → 3'Blocked-NTP



“Two-step-one-pot” synthesis of 3'-blocked nucleotides

Provides scalable, sustainable, economic supply of required ECO Synthesis™ monomers

ECO Synthesis™ Technology: A vision for sustainable RNA synthesis


Accomplished to date

- ✓ Progress on critical TdT performance
- ✓ Proof of Concept for iterative nucleotide addition
- ✓ Concept for enzymatically-derived source of 3'-blocked nucleotides

Next Steps

- ☐ Increase % monomer incorporation to reduce impurities
- ☐ Sustainable & scalable supply of nucleotides
- ☐ Scale-up to process-relevant conditions
- ☐ Achieve gram-scale synthesis of modified RNA

Reimagine enzymes for mRNA manufacturing and RNAi synthesis



Codex® HiCap RNA Polymerase

- ✓ Low dsRNA impurity production
- ✓ High capping efficiency
- ✓ Improved manufacturing economics & efficiency

ECO Synthesis™ Technology

- ✓ New paradigm for synthesis for RNAi therapeutics
- ✓ Sustainable, scalable, economical
- ✓ Proof of concept demonstrated – More to come!