

Quality Control Studies of CF-301 versus *Staphylococcus aureus* ATCC 29213 and *Enterococcus faecalis* ATCC 29212

M. M. TRACZEWSKI ^{1*}, J. Oh² and R. Schuch²

¹ Clinical Microbiology Institute, Wilsonville, OR and ² ContraFect Corp., Yonkers, NY

Maria M. Traczewski BS, MT (ASCP),
Clinical Microbiology Institute
9725 SW Commerce Circle, Ste. A-1
Wilsonville, OR 97070
Tel: 503-682-3232
Email: mtrac@clinmicroinst.com

Abstract

Background: CF-301, a novel, recombinantly-produced bacteriophage-derived lysin (cell wall hydrolase) is the first agent of this class to enter clinical development in the US for the treatment of bacteremia including endocarditis due to *S. aureus*. Previously, a modification of the Clinical and Laboratory Standards Institute (CLSI) broth microdilution method (BMD) was approved for susceptibility testing with CF-301 that uses CAMHB supplemented with 25% horse serum and 0.5 mM DL-Dithiothreitol solution (DTT) to eliminate several phenomenon that resulted in inaccurate and non-reproducible results using the standard method. The studies reported here were performed in order to determine BMD quality control (QC) ranges for CF-301 vs. *S. aureus* ATCC 29213 and *E. faecalis* ATCC 29212 using the new modified broth method. Vancomycin was tested as a control for both studies.

Materials and Methods: CLSI M23-A4¹ (2016) Tier 1 and 2 QC studies were performed using the standard reference method BMD minimum inhibitory concentration (MIC) test (M7-A10², 2015) with the following exceptions: CAMHB was supplemented with 25% horse serum and 0.5 mM DTT. For Tier 1 studies 2 labs tested 5 or 8 replicates a day for 5 days in one lot of supplemented CAMHB. For Tier 2 studies frozen BMD panels containing CF-301 diluted in 3 lots of supplemented CAMHB were supplied to 8 labs and each lab tested 1- 4 replicates per day over a minimum of 3 days.

Results: Table 1. Tier 1 and 2 Quality Control Ranges for CF-301

QC Organism	Tier 1 QC Ranges (µg/mL)	Tier 2 Proposed QC Ranges (µg/mL)
<i>S. aureus</i> ATCC 29213	0.25-1	0.25-2
<i>E. faecalis</i> ATCC 29212	16-64	8-64

Conclusions: The Tier 1 study produced 3 dilution ranges for both organisms. Tier 2 studies using the modified BMD MIC test for CF-301 vs. *S. aureus* ATCC 29213 and *E. faecalis* ATCC 29212 produced reproducible results from the 8 laboratories. These results suggested ranges of 0.25-2 µg/mL for *S. aureus* ATCC 29213 and 8-64 µg/mL for *E. faecalis* ATCC 29212. Vancomycin QC ranges were in control for all tests in both studies. Tier 2 ranges will be submitted to the CLSI for inclusion in the M100³ document.

Introduction

CF-301, a novel, recombinantly-produced bacteriophage-derived lysin (cell wall hydrolase) is the first agent of this class to enter clinical development in the US for the treatment of bacteremia including endocarditis due to *S. aureus*.

Key Features of CF-301 lysin:

- Rapid killing of bacteria by cleaving the structure of the bacterial peptidoglycan cell wall
- Synergy with existing antibiotics
- Rapid clearance of biofilms
- Low propensity to develop resistance
- Fully active against resistant bacteria (while sparing good bacteria)

Acknowledgement

This study was sponsored by a grant from ContraFect Corporation

Materials and Methods

Tier 1 Studies:

- CMI and ContraFect each prepared CLSI frozen broth microdilution MIC test panels containing CF-301 and Vancomycin dilutions
- CF-301 dilutions were prepared in CAMHB supplemented with 25% horse serum and 0.5 mM DL-Dithiothreitol (DTT) solution (CAMHB lysin broth).
- Vancomycin dilutions were prepared in CAMHB.
- CMI tested 25 replicates of each QC strains over 5 days at 5 replicates per day.
- Contrafect tested 40 replicates of each QC over 5 days at 8 replicates per day.

Tier 2 Studies:

- CMI prepared CLSI frozen BMD MIC panels containing CF-301 diluted in 3 different manufacturer's lots of CAMHB lysin broth and vancomycin was included as a control agent in one lot of CAMHB
- CMI shipped MIC panels to 7 outside laboratories for a total of 8 labs doing the testing.
- Each laboratory tested 10 replicates of each QC strain over a 3 to 10 days period (minimum 4 replicates per day and no less than 3 days of testing).
- Labs were instructed to perform repeat tests if any vancomycin results were out of control on day of testing. Vancomycin QC ranges were listed per CLSI M100-S26 document.³
- Results were compiled and analyzed at CMI using CLSI M23-A4 guidelines followed by a second analysis using the RangeFinder Method⁴

Results

Figure 1. *S. aureus* ATCC 29213 CF-301 Broth Microdilution Tier 1 MIC Quality Control Proposed Range 0.25 - 1 µg/mL 100% Included




Figure 2. *E. faecalis* ATCC 29212 CF-301 Broth Microdilution Tier 1 MIC Quality Control Proposed Range 16-64 µg/mL 100% Included




Figure 1. *S. aureus* ATCC 29213 Tier 2 CF-301 in CAMHB+25% Horse Serum and 0.5 mM DTT Proposed Range: 0.25-2 µg/ml = 99.2% Included




Figure 3. *E. faecalis* ATCC 29212 Tier 2 CF-301 in CAMHB+25% Horse Serum and 0.5 mM DTT Proposed Range: 8 - 64 µg/ml = 99.2% Included




Table 1. Tier 1 and Tier 2 Quality Control Ranges for CF-301

QC Organism	Tier 1 QC Ranges (µg/mL)	Tier 2 Proposed QC Ranges (µg/mL)
<i>S. aureus</i> ATCC 29213	0.25-1	0.25-2 µg/ml (99.2%)
<i>E. faecalis</i> ATCC 29212	16-64	8-64 µg/ml (99.2%)

Conclusions

- CF-301 Tier 1 and Tier 2 studies MIC quality control ranges were presented to the CLSI Antimicrobial Susceptibility Testing subcommittee in January 2017.
- 4 dilution ranges for both SA and EF were approved as listed in Table 1 Tier 2.
- The use of 8 labs in Tier 2 vs. 2 labs in Tier 1 slightly increased the variability in MIC ranges as expected.
- Tier 2 study results were within the maximum allowed 4 dilutions acceptable for CLSI MIC QC ranges.
- Vancomycin MIC results were in control for all tests in both studies.

References

1. Development of *In Vitro* Susceptibility Testing Criteria and Quality Control Parameters, 4th Edition, January 2016 Clinical Laboratory Standards Institute, Wayne PA.
2. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically; Approved Standard-Tenth Edition, January 2015, Clinical Laboratory Standards Institute, Wayne, PA.
3. Performance Standards for Antimicrobial Susceptibility Testing M100 26th Edition, Clinical and Laboratory Standards Institute, Wayne PA.
4. Turnidge, J., and G. Bordash. 2007. Statistical methods for establishing quality control ranges for antibacterial agents in Clinical and Laboratory Standards Institute susceptibility testing. *Antimicrob. Agents Chemo.* 51: 2483–2488.