

Melissa J. Karau, M.S.¹, Suzannah M. Schmidt-Malan, M.S.¹, Jayawant Mandrekar, Ph.D.², Dario Lehoux, Ph.D.⁴,
Raymond Schuch, Ph.D.⁴, Cara Cassino, M.D.⁴ Robin Patel, M.D.^{1,3}

¹Division of Clinical Microbiology, ²Division of Biomedical Statistics and Informatics, ³Division of Infectious Diseases, Mayo Clinic, Rochester, MN

⁴Contrafact, Yonkers, NY

Contact Information:

Robin Patel, M.D.

Division of **Clinical Microbiology**

Division of Infectious Diseases

Mayo Clinic, Rochester, MN

Email: patel.robin@mayo.edu

Background

Prosthetic Joint Infection (PJI)^{1,2}

- *Staphylococcus aureus* most common cause of PJI
- Forms biofilms → evades antimicrobials and immune system
- Surgical intervention and long term antibiotic therapy
- Increase in multi-drug resistance

Exebacase

- Phage derived lysin, recombinantly produced as purified cell wall hydrolase enzyme
- Hydrolyzes cell wall of *S. aureus* and associated biofilms
- Rapidly bactericidal, synergistic activity with daptomycin and vancomycin³⁻⁴
- Bone levels ~15% plasma levels after single dose of 10 mg/kg
- Mouse bacteremia and rat osteomyelitis showed decrease in infection with single dose^{3,5}
- Phase 2 clinical trial on bacteremia demonstrated 43% higher clinical response rates *versus* standard therapy alone against MRSA⁶

Treatment

Treatment Groups

- 20mM L-histidine, 5% D-sorbitol (exebacase vehicle)
- Saline (daptomycin vehicle)
- Exebacase 0.098 mg/ml
- Exebacase 0.98 mg/ml
- Exebacase 9.8 mg/ml
- Daptomycin 0.098 mg/ml
- Daptomycin 0.098 mg/ml
- Daptomycin 0.098 mg/ml
- Daptomycin 0.098 mg/ml
- Exebacase and daptomycin each at 0.098 mg/ml


Treatment Durations

- 2 hours
- 4 hours
- 8 hours
- 12 hours

Results

Mean log ₁₀ cfu/K-wire reduction compared to carrier solution								
Hours of Treatment	Daptomycin (DAP)			Exebacase (EXE)			Daptomycin + Exebacase	
	Concentration in mg/ml							
	0.098	0.98	9.8	0.098	0.98	9.8		0.098
2	0.23	0.06	0.46	5.12	4.83	5.68	5.29	
4	0.48	1.21	0.68	3.60	3.20	5.52	3.24	
8	0.39	0.39	0.31	2.24	3.67	5.78	4.34	
12	0.38	0.61	0.81	3.05	4.41	4.76	3.50	

Bold values considered bactericidal



Methods

- MRSA (IDRL-6169) from patient with PJI; minimum inhibitory concentration is 0.5 µg/ml for both daptomycin and exebacase
- Biofilms grown on 5 mm stainless steel threaded Kirschner wires (K-wires) in 40 µl of 10⁶ cfu/ml in tryptic soy broth (TSB) for 10 hours
- Wires removed from growth media → placed in 40 µl of treatment solution
- After specified time, K-wires removed from solution and dipped in sterile saline to remove remnant treatment
- K-wires placed in 0.5 ml sterile saline, vortexed 30 seconds, sonicated (40 kHz, 0.22 W/cm²) 5 minutes, vortexed additional 30 seconds
- Resultant sonicate fluid cultured using 1:10 dilutions in saline, 100 µl of each dilution plated on sheep blood agar → 37°C for 48 hours
- TSB placed into sonicate fluid and also incubated
- Colonies counted or growth from broth only assigned 0.65 log₁₀ cfu/K-wire; negative broth cultures assigned 0.13 log₁₀ cfu/K-wire
- Results reported as log₁₀ cfu/K-wire reduction relative to vehicle alone
- 3-log₁₀ cfu/K-wire reduction considered bactericidal
- Testing performed in triplicate; P-values calculated using Kruskal Wallis

Conclusions

- Exebacase showed rapid activity against MRSA biofilms on orthopedic K-wires
- Exebacase had superior performance to daptomycin
- The local application of exebacase could be a treatment option for the treatment of biofilm related device infections

References

1. Liu C, Chambers HF, Kaplan SL, et al. 2011. Clinical practice guidelines by the Infectious Diseases Society of America for the treatment of methicillin-resistant *Staphylococcus aureus* infections in adults and children. Clin Infect Dis 52:e18-e55.
2. Kavanagh N, Ryan EJ, Widaa A, et al. 2018. Staphylococcal osteomyelitis: Disease progression, treatment challenges, and future directions. Clin Microbiol Rev 31:e00084-17.
3. Schuch R, Lee HM, Schneider BC, Sauve KL, et al. 2013. Combination therapy with lysin CF-301 and antibiotic is superior to antibiotic alone for treating methicillin-resistant *Staphylococcus aureus*-induced murine bacteremia. J Infect Dis 209:1469-1478.
4. Schuch R, Khan BK, Raz A, et al. 2017. Bacteriophage lysin CF-301: A potent anti-staphylococcal biofilm agent. Antimicrob Agents Chemother 61: 02666-16.
5. Karau MJ, Schmidt-Malan SM, Yan Q, Greenwood-Quaintance KE, Mandrekar J, Lehoux D, Schuch R, Cassino C, Patel R. Exebacase in addition to daptomycin is more active than daptomycin or exebacase alone in methicillin-resistant *Staphylococcus aureus* osteomyelitis in rats. ASM Microbe 2019. San Francisco, CA.
6. Fowler VG, Das A, Lipka J, Schuch R, Cassino C. 2019. Exebacase (lysin CF-301) improved clinical responder rates in methicillin-resistant *Staphylococcus aureus* bacteremia and endocarditis compared to standard of care antibiotics alone in a first-in-patient phase 2 study. 29th European Congress of Clinical Microbiology and Infectious Diseases Annual Meeting, Amsterdam, Netherlands.