



BEST IN CLASS SERVICE LIFE

Up to 8 Years*



Inogen®
Rove 6™
System

Inogen® is committed to providing pioneering portable oxygen technology and excellent customer service. Experience best-in-class service life of up to 8 years* with the Inogen® Rove 6™.



DURABILITY

The Inogen® Rove 6™ is tested to withstand normal use for up to 8 years of daily operation and is designed to provide patients the ability to perform their own device maintenance.



FLEET LONGEVITY

Extend the longevity of your concentrator device fleet with up to 8 years of service life, helping to impact the profitability of your oxygen concentrator business.



RELIABILITY

With Inogen's patented Intelligent Delivery Technology, the Rove 6™ is designed to provide the necessary oxygen for most patients during all phases of daily activity and during sleep.

To learn more, please call your local Inogen® representative at
+31 30 7820689 or email us at Info-eu@inogen.net

*When used for less than 8 hours per day.

PM230399 EN_EX_INTL | Rx Only. © 2023 Inogen, Inc.
301 Coromar Drive, Goleta, CA 93117
Inogen® and Rove 6™ are trademarks of Inogen, Inc. The usage of any Inogen, Inc. trademark is strictly forbidden without the prior consent of Inogen, Inc. All other trademarks are trademarks of their respective owners or holders.

USES: The Inogen® Portable Oxygen Concentrator provides a high concentration of supplemental oxygen to patients requiring respiratory therapy on a prescriptive basis. It may be used in home, institution, vehicle, and various mobile environments. **DO NOT USE IF:** This device is not intended to be used in any way other than described in the indications for use. Do not use in parallel or series with other oxygen concentrators or oxygen therapy devices. This device is to be used as an oxygen supplement and is not intended to be life sustaining or life supporting. **ONLY** use this product if the patient is capable of spontaneous breath, able to inhale and exhale without the use of a machine. The conserving, or pulse dose, oxygen delivery technique used by this device is contraindicated in persons whose breathing during normal resting would be unable to trigger the device. Proper device triggering, setup and operation must be confirmed by an experienced clinician or other respiratory professional. Not for pediatric use. Not for use by tracheotomized patients. **WARNINGS:** The device produces enriched oxygen gas, which accelerates combustion. Do not allow smoking or open flames within 2m (6.56ft) of this device while in use. If you feel ill or uncomfortable, or if the concentrator does not signal an oxygen pulse and you are unable to hear and/or feel the oxygen pulse, consult your equipment provider and/or your physician immediately. If you are unable to communicate discomfort, you may require additional monitoring and or a distributed alarm system to convey the information about the discomfort and or the medical urgency to your responsible caregiver to avoid harm. Use only spare parts recommended by the manufacturer to ensure proper function and to avoid the risk of fire and burns. To avoid danger of choking or strangulation hazard, keep cords away from children and pets. **TALK TO YOUR HEALTH CARE PROVIDER:** The oxygen flow setting must be determined and recorded for each patient individually by the prescriber, including the configuration of the device, its parts, and the accessories. It is the responsibility of the patient to periodically reassess the setting(s) of the therapy for effectiveness. The proper placement and positioning of the prongs of the nasal cannula in the nose is critical for oxygen to be delivered.

The Inogen® portable oxygen concentrators are available by Rx only. For safe and proper use of Inogen devices, please refer to <https://provider.inogen.com/en/support/product-materials>