

Non-Invasive Blood Glucose Monitoring in People with Diabetes Using an RF Sensor and Venous Blood Comparator

D. Klyve¹, J. Anderson², K. Currie², C. Ward³, K. Pandya³, V. Somers⁴

¹Department of Mathematics, Central Washington University, Ellensburg, WA, USA, ²Know Labs Inc., Seattle, WA, USA, ³Edge Impulse, Inc., San Jose, CA, USA, ⁴Mayo Clinic, Rochester, MN, USA

KNOW LABS

BACKGROUND & AIM

The high expense and discomfort of current blood glucose (BG) monitoring techniques necessitate the development of a continuous, economical, non-invasive monitor. Despite significant efforts in the development of non-invasive BG monitoring solutions, it remains challenging to deliver an accurate, real-time BG measurement. Our ongoing clinical study assesses the accuracy of the novel Know Labs radiofrequency (RF) sensor for non-invasive BG measurement in people with prediabetes and Type 2 diabetes using venous blood as a comparative reference.

METHODS

- The study employed a novel RF sensor that rapidly sweeps frequencies from 500 MHz to 1500 MHz.
- The sensor scanned 10 participants’ forearms continuously over 21 three-hour Oral Glucose Tolerance Tests (75g).
- Venous blood was collected using a peripheral intravenous catheter (PIVC) every five minutes and analyzed using a blood glucose monitoring test system (StatStrip, Nova Biomedical) as reference values.
- Data were preprocessed using smoothing techniques after which an 80/20 split was performed to create model training and held-out test datasets.
- A light Gradient Boosting Machine Learning Model (lightGBM) was trained on 520 paired observations (RF data and venous BG values), then tested on 130 held-out paired observations.

RESULTS

On the held-out test dataset, BG was estimated with a Mean Absolute Relative Difference (MARD) of $11.1 \pm 2.1\%$ relative to venous blood (Table 1). We observed similar accuracy in normoglycemic ($11.0 \pm 2.7\%$) and hyperglycemic ranges ($11.5 \pm 3.1\%$). A Surveillance Error Grid analysis of model accuracy showed 82.3% of estimations in Risk Grade A and 17.7% in Risk Grade B (Figure 1). No estimations fell in the higher Risk Grades.

Glucose Range (mg/dL)	n	MARD (%)	±15%	±20%
Hypoglycemic (<70)	4	9.5 ± 8.3	75.0 ± 4.2	100.0 ± 0.0
Normoglycemic (70-180)	99	11.0 ± 2.7	75.8 ± 0.8	83.8 ± 0.7
Hyperglycemic (>180)	27	11.5 ± 3.1	66.7 ± 1.8	85.2 ± 1.3
Total	130	11.1 ± 2.1	73.8 ± 0.8	84.6 ± 0.6

Table 1: MARD values and percentages falling within 15% and 20% of the reference value by glycemic status. Error values on the MARD give the 95% *t*-Confidence interval. Error bars on the ±15% and ±20% give the 95% *z*-Confidence interval for proportions.

CONCLUSIONS

These interim results suggest that this ML model applied on data from the Know Labs RF sensor can measure BG non-invasively.

Further data collection and model refinement will continue.




Figure 1: Surveillance Error Grid analysis comparing the 130 ML model estimations in the test dataset to the venous blood reference.

Presenting Author Contact Details: Virend Somers, M.D., Ph.D.
somers.virend@mayo.edu | Cardiovascular Medicine, Mayo Clinic, Rochester, Minnesota 55905 | (507) 255-1144

DK and VS are consultants for and own stock in Know Labs. JA is employed by and has stock options in Know Labs. CW and KP are consultants for Know Labs.