

September 9, 2020

New Head-to-Head Study Demonstrates Equivalency and Accuracy of COVID-19 Rapid Point-of-Care Tests, Despite Different Test Sensitivity Claims

Results show clinical study design may influence performance estimates and published sensitivity claims of SARS-CoV-2 antigen tests

FRANKLIN LAKES, N.J., Sept. 9, 2020 /PRNewswire/ -- BD (Becton, Dickinson and Company) (NYSE: BDX), a leading global medical technology company, today announced the results of a new research study that demonstrated performance of the BD Veritor™ System for Rapid Detection of SARS-CoV-2 Assay was equivalent to that of the Quidel Sofia 2 SARS Antigen Fluorescent Immunoassay (AFI).



In the blinded, head-to-head study, 361 patient samples representing a cross-sectional population were tested for the detection of SARS-CoV-2 using both the BD Veritor™ Plus system and the Quidel Sofia 2 SARS AFI. The two testing platforms had a high degree of agreement on testing results. When there was disagreement between the tests (i.e. one showed positive, one showed negative), the samples were then tested using molecular PCR (polymerase chain reaction), which has higher sensitivity than antigen tests. The PCR analysis of the discrepancies favored the BD Veritor™ Plus system, which had more agreement with PCR analysis than the Quidel Sofia 2 SARS AFI. Of the seven discrepancies, PCR analysis agreed with the BD Veritor™ Plus system on five results and agreed with the Sofia 2 SARS AFI on two results. Overall, the study suggests a high degree of agreement in testing accuracy between the two systems.

"The difference in published sensitivity claims for the Quidel Sofia SARS AFI (96.7%) versus the BD Veritor™ Plus system (84%) was not supported by this study, probably due to the differences in study designs and patient populations in this study versus the Sofia EUA study," said Dr. Charles Cooper, M.D., vice president for Medical and Scientific Affairs at BD and a co-author of the study. "This study was designed to assess performance in samples

from patient populations who commonly seek testing in community-based point-of-care locations where rapid antigen tests would typically be most often used. Other studies enrolled patients in emergency rooms and hospitalized patients which, among other study design factors, have the potential to influence antigen testing positive percent agreement¹. The study also suggests that the number of symptoms in selection criteria may influence positive percent agreement of antigen-based point-of-care testing as well."

Dr. Stephanie N. Taylor, M.D., Professor of Medicine, Section of Infectious Diseases at Louisiana State University and a study co-author said, "This study shows that EUA labeled sensitivity and specificity should not be the sole basis for understanding antigen assay performance. Other factors, such as clinical trial enrollment criteria, study population characteristics, clinical setting, reference method sensitivity, overall prevalence in the trial population and the limit of detection for each antigen test are all important considerations in order to arrive at a fuller understanding of overall performance."

In addition to showing a high degree of agreement between the two tests, the study also supports that the combination of reliable accuracy and speed to result (about 15 minutes) of rapid antigen tests are valuable tools to combat the COVID-19 pandemic, particularly in decentralized health care settings where larger-scale, molecular-based tests are not practical or would result in significantly longer turnaround times.

As part of its emergency use authorization (EUA) for the SARS-CoV-2 assay for the BD Veritor™ Plus system, BD completed one of the geographically broadest U.S. clinical trials of any antigen test on the market today, collecting and analyzing 226 samples from 21 diverse areas across the United States. This effort was intended to reflect performance across a range of different community-based clinical settings by enrolling patients reflective of those typically seeking rapid point-of-care diagnostic testing. The BD Veritor™ Plus system was studied in care settings and patient populations most likely to rely on rapid point-of-care tests, including ambulatory care sites, outpatient clinics, drive-through and walk-up testing sites, and skilled nursing facilities.

The research study authors have submitted the paper for peer-review consideration to an infectious disease journal. The preprint article titled "Clinical evaluation of BD Veritor™ SARS-CoV-2 point-of-care test performance compared to PCR-based testing and versus the Sofia SARS Antigen point-of-care test" is currently available on [medRxiv](https://www.medrxiv.org/content/10.1101/2020.07.27.20161111v1).

1. Positive percent agreements shown in the EUA instructions for use are estimates of sensitivity.

About the BD Veritor™ System for Rapid Detection of SARS-CoV-2 Assay

The BD Veritor™ Plus System for Rapid Detection of SARS-CoV-2 Assay has not been cleared or approved by FDA. The test has been authorized by FDA under an EUA for use by authorized laboratories. The test has been authorized only for the detection of proteins from SARS-CoV-2, not for any other viruses or pathogens; and, the test is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b) (1) of the Act, 21 U.S.C. § 360bbb-3(b)(1), unless the authorization is terminated or revoked sooner.

About BD

BD is one of the largest global medical technology companies in the world and is advancing the world of health by improving medical discovery, diagnostics and the delivery of care. The

company supports the heroes on the frontlines of health care by developing innovative technology, services and solutions that help advance both clinical therapy for patients and clinical process for health care providers. BD and its 65,000 employees have a passion and commitment to help enhance the safety and efficiency of clinicians' care delivery process, enable laboratory scientists to accurately detect disease and advance researchers' capabilities to develop the next generation of diagnostics and therapeutics. BD has a presence in virtually every country and partners with organizations around the world to address some of the most challenging global health issues. By working in close collaboration with customers, BD can help enhance outcomes, lower costs, increase efficiencies, improve safety and expand access to health care. For more information on BD, please visit bd.com.

Contacts:

Troy Kirkpatrick
BD Public Relations
858.617.2361
troy.kirkpatrick@bd.com

Kristen Stewart
BD Investor Relations
201.847.5378
kristen.stewart@bd.com



View original content to download multimedia:<http://www.prnewswire.com/news-releases/new-head-to-head-study-demonstrates-equivalency-and-accuracy-of-covid-19-rapid-point-of-care-tests-despite-different-test-sensitivity-claims-301126072.html>

SOURCE BD (Becton, Dickinson and Company)