


# FREEDOM60<sup>®</sup>

Part of the Freedom Integrated Syringe Infusion System



## 123

### Easy to use & learn

Only 2 controls.  
Just wind, load & go!



### Safe constant pressure

Highly accurate constant pressure system, meeting ISO accuracy standards for safe drug delivery



### DynEQ™ (Dynamic Equilibrium)

If the pressure in your sites increases, the system will immediately respond by automatically slowing down.



### Quality assured

Made in the USA with the highest level of quality control for defect-free products.



### Advanced technology

All parts are designed to fit & work together within the syringe driver. The Negator motor is the 'heart' of the FREEDOM60<sup>®</sup>.



### Portable

No batteries or electricity needed. Comes with a travel pouch to promote patient independence.



### Fits 50 ml BD<sup>®</sup> syringe

Easy to handle and purge of air. Can be filled from 1 ml to 50 ml.



### Home use

No sensitive electric parts or preventive maintenance, made to withstand frequent home use.



### Spring-loaded syringe holder

Locks the syringe into the driver to maintain proper syringe orientation, prevents inadvertent removal during administration.



### Large wind knob

Ergonomic, easy to grasp & turn, disengages when the FREEDOM60<sup>®</sup> is ON, cannot be over-wound.



### Recessed on/off switch

No unintentional on or off.



### Syringe shield

Easy to check infusion progress, protects the syringe.

**NOTE:** Use ONLY KORU Freedom System accessories with the FREEDOM60<sup>®</sup> Syringe Driver. Use of any other accessory may cause over delivery of fluids, unknown flow rates or additional site complications such as pain, swelling and redness.

## Related Freedom Integrated System Products:



### FreedomEdge<sup>®</sup> Syringe Infusion System

For smaller infusion doses, BD<sup>®</sup> 20 & 30 ml syringes



### HiGH-Flo Subcutaneous Safety Needle Sets™

24G & 26G, 4 -14 mm lengths, 1-6 leg sets



### Precision Flow Rate Tubing™

For flow rates from F0.5 - F2400



### 24" Extension Tubing Set (FEXT)

For those who need longer tubing length



KORU Medical Systems | 845-469-2042 | korumedical.com | @korumedical

The FREEDOM60<sup>®</sup> and FreedomEdge<sup>®</sup> Syringe Infusion Systems, Precision Flow Rate Tubing™ and HiGH-Flo Subcutaneous Safety Needle Sets™ are trademarks of KORU Medical Systems and are compliant with Medical Device Directive 93/42/EEC. KORU Medical Systems is ISO 13485 certified. ©2020 KORU Medical Systems; All Rights Reserved. F60\_FL\_Benefits\_03-10005vD