


# HighFlo™ Subq Safety Needles

Part of the Freedom Integrated Syringe Infusion System




HighFlo



other brand


## Backcut needle tip design

Minimizes damage, scarring & local site reactions.



## Flexible flying hinge

Wings have great flexibility & allow each needle to stay in place without tenting.



## Safety closure

Protects against accidental needle stick injuries.




## Custom luer

Designed to provide even flow to every site.



## Gauge options

Designed for subq flow rates with options in 26 & 24 gauge.



## Length options

Needles are available in 4, 6, 9, 12, & 14 mm lengths.



## Tubing length

Efficient & practical length for less tangles & faster flow rates.  
20" in HighFlo 26G sets;  
24" in HighFlo 24G.



## Y-connector


Combines sets for more than 6 sites (up to 8), or to enable infusing with different length needles.




## Needle site dressing

3M Tegaderm™ comfort adhesive technology, strong hold & gentle removal, included with every needle set.

## 26G set options & residual volume



## 24G set options & residual volume



## Related Freedom Integrated System Products:



### FREEDOM60® Syringe Infusion System

For larger infusion doses, BD® 50 ml syringe



### FreedomEdge® Syringe Infusion System

For smaller infusion doses, BD® 20 & 30 ml syringes



### Precision Flow Rate Tubing™

For flow rates from F0.5 - F2400



### 24" Extension Tubing Set (FEXT)

For those who need longer tubing length



KORU Medical Systems | 845-469-2042 | korumedical.com | @korumedical

The FREEDOM60® and FreedomEdge® Syringe Infusion Systems, Precision Flow Rate Tubing™ and HighFlo Subcutaneous Safety Needle Sets™ are trademarks of KORU Medical Systems and are compliant with Medical Device Directive 93/42/EEC. KORU Medical Systems is ISO 13485 certified. ©2020 KORU Medical Systems; All Rights Reserved. HFL\_FL\_Benefits\_03-10003vF