



Stellar KLH™
ELISA Test Kit

Stellar
BIOTECHNOLOGIES


Keyhole Limpet Hemocyanin (KLH)

ELISA Immunoassay

Rapid
Specific
Sensitive
Quantitative

Stellar Biotechnologies Anti-KLH ELISA Test Kits. Highly sensitive Enzyme-Linked Immunosorbent Assays designed for the rapid, quantitative determination of anti-KLH antibodies from animal serum or plasma.

WWW.STELLARBIOTECH.COM



The Only ELISA Made with **Stellar KLH**™

Keyhole Limpet Hemocyanin (KLH) is a large, oxygen-carrying glycoprotein from the marine mollusk *Megathura crenulata*. KLH is well known as a safe, potent stimulator of humoral and cellular immune responses. It plays a vital role in research and clinical studies as an antigen for assessing immune function such as in monitoring the immunosuppressive effects of drug candidates.



Stellar KLH™ is made by Stellar Biotechnologies, the world leader in high quality and sustainable KLH technology. Our KLH manufacturing and experience are unmatched. We developed industry-leading methods that ensure quality, consistent protein with optimum immunogenicity for your research needs, while protecting the natural source of this important biomedical material.

UNIQUE BENEFITS

The quality of KLH used for solid phase (microtiter well) immobilization can affect test results. Only Stellar brand assay plates are pre-coated with **Stellar KLH**™ protein which provides wide dynamic range, low assay background noise and quantitative measurement.

- Quantitative ELISA
- Wide Dynamic Range
- Low Assay Background
- Reproducible Linearity
- Validated with Controls
- Consistent and Accurate Results



Broad Suite of Anti-KLH ELISA Kits

Stellar KLH™ ELISA test kits are available to measure either IgG or IgM, in non-human primates, mice and rat models. Contact Stellar Biotechnologies for availability and pricing.




Stellar Biotechnologies, Inc. is a world leader in the development and manufacture of Keyhole Limpet Hemocyanin (KLH).

ORDERING INFORMATION (CALL FOR PRICING)

Product Code	Description
ELI-01G	Mouse Anti-KLH IgG ELISA Test Kit
ELI-01M	Mouse Anti-KLH IgM ELISA Test Kit
ELI-02G	Rat Anti-KLH IgG ELISA Test Kit
ELI-02M	Rat Anti-KLH IgM ELISA Test Kit
ELI-03G	Monkey Anti-KLH IgG ELISA Test Kit
ELI-03M	Monkey Anti-KLH IgM ELISA Test Kit


RELATED PRODUCTS (CALL FOR PRICING)

Stellar KLH™ protein, for quality and consistent results.  Stellar KLH™ Protein

Description

Native Stellar KLH	High Molecular Weight KLH (+8MDa) in DPBS
Subunit Stellar KLH	Subunit formulations (+440 – 800 KDa) in WFI

Procedure



Stellar

BIOTECHNOLOGIES

Stellar Biotechnologies, Inc.
332 East Scott Street, Port Hueneme, CA 93041 USA
Toll Free: (855) KLH-7555 / +1 (805) 488-2800
Email: KLHinfo@stellarbiotech.com
www.stellarbiotech.com

Stellar Biotechnologies, Inc. is an immuno-pharmaceutical company and the world leader in development and manufacture of high-quality, sustainable Keyhole Limpet Hemocyanin (KLH).

We revolutionized the KLH industry with production know-how and environmentally sound methods which turned this vital, natural protein into a renewable, quality one. Stellar Biotechnologies is committed to providing excellence in targeted markets ranging from KLH-enhanced products to therapeutic vaccines and immunodiagnostics.



Stellar KLH™
ELISA Test Kit

THREE WAYS TO ORDER:



Email

KLHInfo@stellarbiotech.com



Fax

+1 (805) 488-1278



Phone

U.S. Toll Free (855) KLH-7555
International +1 (805) 488-2800