B-004

AACR 10/29/17

Enhanced Efficacy of Immune Modulators with Albumin Binding Domains (ABD)


Haomin Huang, Keneneshia Haenssen, Anil Bhate, Supriya Sanglikar, John Baradei, Shan Liu, Senthil Kumar, Zihao Cui, Richard Hampton, Robert Kramer and John Cini Sonnet Biotherapeutics, 1 Duncan Drive, Cranbury, NJ 08512, USA


BioTherapeutics


INTRODUCTION


Small (<100Kd) recombinant therapeutic proteins penetrate better into solid tumors but exhibit shorter circulation half-life (mins/hour vs. 21 days) when compared to IgG's. Methods to improve serum half-life for these small MW therapeutic proteins include conjugation to polyethylene glycol (ie., PEG: eg., Peg-Intron, NKT-214) or formulations on the surface of, or within, nanoparticles. Albumin has also been successfully exploited as a carrier for small therapeutic proteins. Examples include the genetic fusion of Interferon to albumin (Albuferon) and the use of a multivalent albumin binding antibody fragment that has the capacity to bind serum albumin and also neutralize TNFa (Ozoralizumab). A potential advantage to the use of albumin as a drug carrier is that albumin is known to accumulate in inflamed and angiogenic tissues such as the tumor microenvironment, and may thus provide a means of passive targeting to the tumor.


Sonnet BioTherapeutics, using a XOMA phage library, has developed a fully human single chain variable antibody (scFv) construct with demonstrated high binding affinity to serum albumin across species (mouse, human & cynomolgus monkey) and retains the benefits of FcRn mediated recycling of albumin for extending serum half-life. This scFv albumin binding domain (ABD) has been utilized to create several different genetic fusion constructs with various small molecular weight therapeutic proteins (e.g., recombinant interleukins such as IL12, and scFv's against relevant immune-oncology targets, ie., TGFb). In mouse models, these ABD constructs have extended serum half-lives, improved tumor accumulation and enhanced efficacy compared to their respective naked recombinant therapeutic protein. This poster summarizes the discovery of our scFv-ABD and preliminary studies with four different therapeutic scFv-ABD constructs demonstrating the broad utility of this approach.



















Sonnet Tri-Specific Construct


Enhanced Biologic Activity:

IL12: ↑ IL15 alpha receptor, IFN, NK/Tcells, TH1 & ↓ Treg

IL15: ↑ IL12 beta 1 receptor, NK cells & ↓ CD8 memory loss










- . A novel 26 Kd scFv antibody (ABD) that binds mouse, cynomolgus monkey and human serum albumin has been identified, and is being developed, as a platform for enhancing the pharmacokinetics and tumor targeting of small molecular weight therapeutic proteins (Figures 1-3).
- We demonstrate increased retention of an anti-TGF β -scFv-ABD construct in 4T1 tumors in mice (Figure 4) and up to a 10-fold enhancement of the serum half-life in mice of recombinant IL12 or IL15 when fused to the scFv-ABD construct (Figure 5).
- 3. Dose dependent increases in anti-tumor activity were demonstrated with either IL15-ABD or IL12-ABD constructs (Figures 6-8), producing a corresponding increase in tumor infiltrating NK cells for IL15-ABD (Figures 6 & 7). FACs studies on cell infiltrates following IL12-ABD treatment are planned.
- 1. Molar equivalent doses of IL12-ABD were significantly more efficacious (3-30 fold) than naked recombinant IL12 vs. the B16F10 mouse melanoma model (Figure 8 & 10) and produced a corresponding increase in immune response as reflected by spleen weight and serum IFN, which was transient and had no effect on mouse body weight (Figure 9).
- 5. A first-in-class dual targeting tri-specific cytokine construct (IL15-ABD-IL12) demonstrated slightly enhanced in vitro cellular activity compared to either recombinant IL12 or IL15 (Figure 13), that translated to superior in vivo activity compared to IL12-ABD vs. the B16F10 mouse melanoma model (Figure 14).

Sonnet wishes to acknowledge Drs. Olesia Buiakova & Yan Xu (InVivoTek, NJ) for their excellent work involving all animal studies,



www.sonnetbio.com


info@sonnetbio.com


ABD Bio-Activity Selection Criteria to Ensure Increase in PK ☑ Binding to Human & Cyno Serum Albumin kD ~20 - 60 nM \checkmark Binding to Mouse Serum Albumin $kD \sim 10 - 30 \, nM$ ✓ Binding at low pH 5.5 & neutral pH7.2 for human, mouse albumin ✓ No competition with FcRn binding to serum albumin. In-vitro studies with ABD showed cell recycling with the FcRn hence should extend pK in-vivo (see Fig below) Binds Albumin in serum by immune-precipitation ✓ Human serum stability at 37C for 7-14 Days

Immunogenicity tested by Antitope computer analysis XOMA (Library 47) full human scFv phage library was used to screen for a scFv by solid & solution panning against HSA, MSA and h/mTGFβ 1, 2 & 3. ABD selection was based on meeting all of the above criteria. Anti Albumin Binding Domain (ABD) A10M3 is a scFv: VH3/VL1, Tm ~60C, MW 26KD, pl 7.1

Anti TGFβ selection based on blocking TGF inducing of: Tregs proliferation, TGF→SMAD2 PO4, EMT & Migration Anti TGFβ (4D9M) is a scFv: VH3/VL1, Tm ~62C, MW 26KD, pl 5.3 Figure 2

Evaluation of anti-TGF β 4D9M scFv & anti-TGF β 4D9M-ABD Tumor Accumulation In vivo Proof of Concept





ABD accumulates in the tumors 24 hrs+ and without ABD, the scFv enters tumor but diffuses out after 4 hrs.